

Haemorrhoid Multi-Band Ligator

Haemoband-PLUS

Automatic Reloading Technology

Single-Use

Product No. H1002

From:
Haemoband
SURGICAL

The Innovative design includes:

GREEN MOULDED BANDS – These have been developed from our unique composite material specifically designed to ensure consistent performance as well as an extended shelf life. The individually moulded bands roll off more easily, eliminating twisting and ensuring a more reliable release and secure fit onto the tissue.

TRANSPARENT SUCTION NOZZLE – The clarity of the see-through tip maximises visibility throughout the procedure, positioned at the tip of a slim 120mm insertion probe.

PRE-LOADED WITH 4 BANDS – The device is clean and ready-to-use immediately from the packaging, without the need to load any new bands. The device automatically reloads an additional 3 bands for ligation of up to three haemorrhoids, with the fourth band as spare.

MULTI-ACTION HANDLE – The innovative handle controls the suction, fires and reloads the bands. A new stepped line guides the user through each stage of the application.





Gold Standard Treatment

There are a wide range of options available in the treatment of Haemorrhoids; depending on the severity of the condition these range from simple dietary advice to aggressive surgical intervention. Minimally invasive, outpatient procedures are obviously preferred by patients.

Rubber Band Ligation using the Haemoband-PLUS device is highly regarded by patients and clinicians alike, as an efficient and safe treatment for symptomatic 1st -3rd degree haemorrhoids.

The procedure can be performed on an outpatient basis and is proven to be straightforward, easy and economical.



Key Benefits:

Clinically Efficient

Haemoband-PLUS is clinically very efficient. Its multi-action handle means you only need one hand to control suction, release, and reload bands. This leaves the other hand free to hold a proctoscope and allows the user to work alone. The design of Haemoband-PLUS also allows for much more accurate applications, reducing the need for additional procedures.

Easier & Faster

The gun is pre-loaded, already clean and can be used single-handed. The time taken for a typical ligation procedure is reduced from 20 minutes to between 2 and 4 minutes for the banding of up to 3 haemorrhoids. This not only improves patient comfort considerably, but its disposability negates the risk of cross-contamination.

Ready-To-Use

The Haemoband-PLUS gun is ready to use straight from the airtight, clean bag. There are no additional parts, and the only extra equipment you need is a proctoscope and standard hospital suction.

Cost Effective

The Haemoband-PLUS is one of the most cost effective multi-ligators available. It is a Single-Use device and therefore eliminates the additional costs associated with multi-use devices.

How To Use:



Step 1:

Attach suction using the universal applicator at the bottom of the handle.



Step 2:

Press the trigger to the first marked line to induce suction - this will draw the desired tissue into the nozzle.



Step 3:

When tissue is secured, press the trigger to its final position - this will release the band onto the targeted haemorrhoid.



Step 4:

Release! When the trigger is completely released a new band will be automatically reloaded and ready for the next band application.



Haemoband Surgical Ltd.
The Mount Conference Centre
2 Woodstock Link
Belfast
BT6 8DD
Northern Ireland, United Kingdom.

Tel: +44 (0)28 90737283
Email: info@haemobandsurgical.com